

Hemostasis Optimization Program

Continuous Oozing

Will not stop with compression/ simple packing.

The solution for this bleeding is more time consuming than it is difficult.¹

SURGICEL™

Absorbable Hemostats



Problematic

Even though the bleeding is accessible, it could be trouble. It is more than routine and likely to be **resistant to conventional means** and requires immediate attention causing disruption to the normal progression of surgery.¹

EVARREST™

Fibrin Sealant Patch



Difficult to access

Bleeding that occurs **in tight and irregular spaces** and you cannot see the exact source of the bleed. You are concerned accessing a tight space will cause more harm.¹

SURGIFLO™

Hemostatic Matrix



Potential re-bleeding risk

Bleeding may be addressed intraoperatively, but could later **develop into more serious complications**, especially in high-risk patients.¹

VISTASEAL™

Fibrin Sealant (Human)



J&J MedTech

Please see Important Safety Information attached to back.

ETHICON

**Essential Product Information and Important
Safety Information is adhered to the card.
Please see accompanying full Prescribing
Information for complete details.**

The bleeding situations identified reflect customer insights/market research on optimal adjunctive hemostat utilization. The product solutions should only be used in accordance with their instructions for use. Product recommendations should not supplant medical judgment. Surgeon preference, experience, and patient needs may dictate alternate technique. Review all relevant precautions, especially the indications, contraindications, warnings, and information for use. Please see package inserts for Full Prescribing Information. The visual does not reflect any sequential order in use.

References: 1. Hemostasis Optimization Program and Next Steps Meeting, dated 06/30/2014.

SURGICEL™ Powder Absorbable Hemostat Essential Product Information
INDICATIONS

SURGICEL™ Powder (oxidized regenerated cellulose) is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial hemorrhage when ligation or other conventional methods of control are impractical or ineffective. SURGICEL Powder can also be applied in laparoscopic or other endoscopic procedures when used with the SURGICEL™ Endoscopic Applicator. The SURGICEL Endoscopic Applicator is intended for use in delivering SURGICEL Powder absorbable hemostat to bleeding surgical sites through a 5 mm or larger trocar.

CONTRAINDICATIONS

- Do not inject or place SURGICEL Powder into an open blood vessel. Do not use to treat bleeding from large defects in arteries or veins.
- SURGICEL Powder should not be used for implantation in bone defects, such as fractures, since there is a possibility of interference with callus formation and a theoretical chance of cyst formation.
- When SURGICEL Powder is used to help achieve hemostasis in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, or the optic nerve and chiasm, it must always be removed after hemostasis is achieved since it will swell and could exert unwanted pressure.
- SURGICEL Powder should not be used to control hemorrhage from large arteries or veins.
- SURGICEL Powder should not be used on non-hemorrhagic serous oozing surfaces, since body fluids other than whole blood, such as serum, do not react with SURGICEL Powder to produce satisfactory hemostatic effect.
- SURGICEL Powder is an absorbable hemostat, and should not be used as an adhesion prevention product.
- The SURGICEL Powder and the SURGICEL Endoscopic Applicator devices were not designed for intraluminal procedures.

WARNINGS

- SURGICEL Powder is not intended for use on dry (non-bleeding) surfaces or for prevention of bleeding.
- SURGICEL Powder is not intended as a substitute for careful surgery and the proper use of sutures and ligatures.
- Closing with SURGICEL Powder in a contaminated wound without drainage may lead to complications and should be avoided.
- The hemostatic effect of SURGICEL Powder is greater when it is applied dry; therefore, it should not be moistened with water or saline prior to application.
- SURGICEL Powder should not be impregnated with anti-infective agents or with other materials such as buffering or hemostatic substance. Its hemostatic effect is not enhanced by the addition of thrombin, the activity of which is destroyed by the low pH of the product.
- Although SURGICEL Powder may be left in situ when necessary, it is recommended to remove excess powder with irrigation and aspiration once hemostasis is achieved, without disturbing the clot.

- Dislodgement of SURGICEL Powder could possibly occur by intraoperative manipulation, lavage, exaggerated respiration, etc. With other SURGICEL products there have been reports that in procedures such as lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe, when the product was left in the patient after closure it migrated from the site of application into foramina in bone around the spinal cord, resulting in paralysis and, in one case, the product migrated into the left orbit of the eye, causing blindness. While these reports cannot be confirmed to be related to SURGICEL products, special care must be taken by physicians, regardless of the type of surgical procedure. Consider removing SURGICEL Powder in these applications (procedures) after hemostasis is achieved.
- SURGICEL Powder is dry and there may be difficulties in precise delivery under certain circumstances. Unintentional device placement may result in powder scattering and device migration that may increase the risk of adhesion formation. In preclinical in vivo animal studies it was demonstrated that SURGICEL Powder does not increase the incidence of remote adhesions in laparoscopic procedures.
- Although SURGICEL Powder is bactericidal against a wide range of pathogenic microorganisms, it is not intended as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or to prevent postoperative infections.
- To prevent clogging with the SURGICEL Endoscopic Applicator Tip, do not touch the tip to wet surface. Be careful to avoid damaging tissue with the rigid tip.
- Do not attempt to trim the applicator tip. Replace the tip if it becomes clogged.

PRECAUTIONS

- SURGICEL Powder should not be used in conjunction with autologous blood salvage circuits, because its fragments may pass through the transfusion filters of blood-scavenging systems.
- Use minimal amount of SURGICEL Powder required to achieve hemostasis, and remove excess powder in the area of drains to prevent clogging.
- Use only as much SURGICEL Powder (oxidized regenerated cellulose) as is necessary and apply only where needed for hemostasis. Remove any excess before surgical closure in order to facilitate absorption and to minimize the possibility of foreign body reaction, such as encapsulation of the product, which may mimic artifacts on radiographic images, resulting in diagnostic errors and possible reoperation.
- In urological procedures, minimal amounts of SURGICEL Powder should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product.
- Since absorption of SURGICEL Powder could be prevented in chemically cauterized areas, its use should not be preceded by application of silver nitrate or any other escharotic chemicals.
- If SURGICEL Powder is used temporarily to line the cavity of open wounds, it should be removed by irrigation with sterile water or saline solution after bleeding has stopped.
- Precautions should be taken in otorhinolaryngologic surgery to ensure that none of the material is aspirated by the patient (e.g., controlling hemorrhage after tonsillectomy and controlling epistaxis).

- The applicator tip provided on the SURGICEL Powder device is not intended for laparoscopic or other endoscopic use. If laparoscopic or other endoscopic use is desired, remove the existing applicator tip from the SURGICEL Powder device, and replace with the SURGICEL Endoscopic Applicator tip (supplied separately). In laparoscopic or other endoscopic procedures, SURGICEL Powder should only be applied using the SURGICEL Endoscopic Applicator. Consult the SURGICEL Endoscopic Applicator Instructions for Use (IFU) for proper assembly and directions for use with the SURGICEL Powder device.
- The SURGICEL Endoscopic Applicator is supplied with a flexible inner tip inside a rigid cannula. The rigid cannula cannot be used independently.
- The SURGICEL Endoscopic Applicator should only be used by persons having adequate training and familiarity with endoscopic techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any endoscopic procedure.
- To prevent inadvertent device spillage, or unintended contact with tissue, organs, or blood, maintain visualization of the SURGICEL Endoscopic Applicator tip at all times.
- Do not compress or excessively bend the flexible inner tip of the SURGICEL Endoscopic Applicator which could obstruct the application of the powder. It is possible that the powder accumulated in the applicator could disperse beyond the target bleeding site upon compression of the bellows, which may require additional irrigation and aspiration.

ADVERSE EVENTS

- Paralysis and nerve damage have been reported when other SURGICEL products were used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm.
- Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when other SURGICEL products were placed in the anterior cranial fossa (see WARNINGS and PRECAUTIONS).
- Foreign body reactions have been reported with other products from the SURGICEL Family of Absorbable Hemostats. Burning has been reported when other SURGICEL products were applied after nasal polyp removal. Headache, burning, stinging, and sneezing in epistaxis and other rhinological procedures, and stinging when SURGICEL product was applied on surface wounds (varicose ulcerations, dermabrasions, and donor sites) have also been reported.
- Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy have been reported.

For more information and technical questions, call 1-800-795-0012. For complete information including indications, contraindications, warnings, precautions, adverse reactions, and directions for use, consult the product package insert.

US_ETH_BIOS_109671

SURGICEL™ Essential Product Information
INDICATIONS

SURGICEL™ Absorbable Hemostat (oxidized regenerated cellulose) is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial hemorrhage when ligation or other conventional methods of control are impractical or ineffective. SURGICEL™ ORIGINAL, SURGICEL™ FIBRILLAR™, SURGICEL NU-KNIT™, and SURGICEL SNoW™ Absorbable Hemostats can be cut to size for use in endoscopic procedures.

CONTRAINDICATIONS

- Although packing or wadding sometimes is medically necessary, SURGICEL Absorbable Hemostat should not be used in this manner, unless it is to be removed after hemostasis is achieved (See WARNINGS and PRECAUTIONS).
- SURGICEL Absorbable Hemostat should not be used for implantation in bone defects, such as fractures, since there is a possibility of interference with callus formation and a theoretical chance of cyst formation.
- When SURGICEL Absorbable Hemostat is used to help achieve hemostasis in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, or the optic nerve and chiasm, it must always be removed after hemostasis is achieved since it will swell and could exert unwanted pressure.
- SURGICEL Absorbable Hemostat should not be used to control hemorrhage from large arteries.
- SURGICEL Absorbable Hemostat should not be used on non-hemorrhagic serous oozing surfaces, since body fluids other than whole blood, such as serum, do not react with SURGICEL Absorbable Hemostat to produce satisfactory hemostatic effect.
- SURGICEL Absorbable Hemostat is an absorbable hemostat and should not be used as an adhesion prevention product.

WARNINGS

- SURGICEL Absorbable Hemostat is not intended as a substitute for careful surgery and the proper use of sutures and ligatures.
- Closing SURGICEL Absorbable Hemostat in a contaminated wound may lead to complications and should be avoided.
- The hemostatic effect of SURGICEL Absorbable Hemostat is greater when it is applied dry; therefore it should not be moistened with water or saline.
- SURGICEL Absorbable Hemostat should not be impregnated with anti-infective agents or with other materials such as buffering or hemostatic

- substances. Its hemostatic effect is not enhanced by the addition of thrombin, the activity of which is destroyed by the low pH of the product.
- Although SURGICEL Absorbable Hemostat may be left in situ when necessary, it is advisable to remove it once hemostasis is achieved. It must always be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm, and in proximity to tubular structures that could become constricted by swelling, regardless of the type of surgical procedure because SURGICEL Hemostat, by swelling, may exert pressure resulting in paralysis and/or nerve damage. Dislodgement of SURGICEL Absorbable Hemostat could possibly occur by means such as repacking, further intraoperative manipulation, lavage, exaggerated respiration, etc. There have been reports that in procedures such as lobectomy, laminectomy and repair of a frontal skull fracture and lacerated lobe that SURGICEL Absorbable Hemostat, when left in the patient after closure, migrated from the site of application into foramina in bone around the spinal cord resulting in paralysis and, in another case, the left orbit of the eye, causing blindness. While these reports cannot be confirmed, special care must be taken by physicians, regardless of the type of surgical procedure, to consider the advisability of removing SURGICEL Absorbable Hemostat after hemostasis is achieved.
 - Although SURGICEL Absorbable Hemostat is bactericidal against a wide range of pathogenic microorganisms, it is not intended as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or prevent post-operative infections.

PRECAUTIONS

- Use only as much SURGICEL Absorbable Hemostat as is necessary for hemostasis, holding it firmly in place until bleeding stops. Remove any excess before surgical closure in order to facilitate absorption and minimize the possibility of foreign body reaction, such as encapsulation of the product, which may mimic artifacts on radiographic images, resulting in diagnostic errors and possible reoperation.
- In urological procedures, minimal amounts of SURGICEL Absorbable Hemostat should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product.
- Since absorption of SURGICEL Absorbable Hemostat could be prevented in chemically cauterized areas, its use should not be preceded by application of silver nitrate or any other escharotic chemicals.

- If SURGICEL Absorbable Hemostat is used temporarily to line the cavity of large open wounds, it should be placed so as not to overlap the skin edges. It should also be removed from open wounds by forceps or by irrigation with sterile water or saline solution after bleeding has stopped.
- Precautions should be taken in otorhinolaryngologic surgery to ensure that none of the material is aspirated by the patient. (Examples: controlling hemorrhage after tonsillectomy and controlling epistaxis.)
- Care should be taken not to apply SURGICEL Absorbable Hemostat too tightly when it is used as a wrap during vascular surgery (see Adverse Reactions section of the complete product package insert).

ADVERSE EVENTS

- Paralysis and nerve damage have been reported when SURGICEL Absorbable Hemostat was used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm.
- Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when SURGICEL Absorbable Hemostat was placed in the anterior cranial fossa
- "Encapsulation" of fluid and foreign body reactions have been reported.
- Burning has been reported when SURGICEL products were applied after nasal polyp removal. Headache, burning, stinging, and sneezing in epistaxis and other rhinological procedures, and stinging when SURGICEL product was applied on surface wounds (varicose ulcerations, dermabrasions, and donor sites) have also been reported.
- There have been reports of stenotic effect when SURGICEL Absorbable Hemostat has been applied as a wrap during vascular surgery.
- Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy have been reported.

For more information, please consult your doctor or for product quality and technical questions, call 1-800-795-0012. For complete product information including indications, contraindications, warnings, precautions, and adverse reactions, please reference the individual product package inserts.

US_ETH_BIOS_109987

SURGIFLO™ Hemostatic Matrix Kit Essential Product Information (Made from Absorbable Gelatin Sponge, USP) with Thrombin
DESCRIPTION
SURGIFLO™ with Thrombin (SURGIFLO™ Hemostatic Matrix Kit) is intended for hemostatic use by applying to a bleeding surface.

ACTIONS
When used in appropriate amounts SURGIFLO™ is absorbed completely within 4 to 6 weeks.

INTENDED USE/INDICATIONS
SURGIFLO™, mixed with thrombin solution, is indicated in surgical procedures (other than ophthalmic) as an adjunct to hemostasis when control of bleeding by ligature or other conventional methods is ineffective or impractical.

CONTRAINDICATIONS

- Do not use SURGIFLO™ in intravascular compartments because of the risk of embolization.
- Do not use SURGIFLO™ in patients with known allergies to porcine gelatin.
- Do not use SURGIFLO™ in closure of skin incisions because it may interfere with the healing of skin edges. This interference is due to mechanical interposition of gelatin and is not secondary to intrinsic interference with wound healing.

WARNINGS

- SURGIFLO™ Hemostatic Matrix is not intended for prevention of bleeding. The use of SURGIFLO™ Hemostatic Matrix for mechanical support has not been studied.
- SURGIFLO™ should not be used in the presence of infection and should be used with caution in contaminated areas of the body.
- SURGIFLO™ should not be used in instances of pumping arterial hemorrhage. SURGIFLO™ will not act as a tampon or plug in a bleeding site.
- SURGIFLO™ should be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm because it may swell resulting in nerve damage.
- Excess SURGIFLO™ should be removed once hemostasis has been achieved.
- The safety and effectiveness of SURGIFLO™ for use in ophthalmic procedures has not been established.
- SURGIFLO™ should not be used for controlling post-partum intrauterine bleeding or menorrhagia.
- The safety and effectiveness of SURGIFLO™ has not been established in children and pregnant women.
- The blue flexible applicator tip should not be trimmed to avoid exposing internal guidewire.
- The white straight applicator tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip.

PRECAUTIONS

- Safe and effective use of SURGIFOAM™ Sponge has been reported in a published neurologic retrospective study involving 1700 cases in Europe. Safe and effective use in neurosurgery has not been proven through

- randomized, controlled clinical studies in the United States.
- SURGIFLO™ Hemostatic Matrix is for single use only. Do not resterilize. If the product is reused the performance of the product may deteriorate, and cross contamination may occur which may lead to infection.
- SURGIFLO™ Hemostatic Matrix is supplied as a sterile product. Unused open SURGIFLO™ Hemostatic Matrix should be discarded. Do not use SURGIFLO™ Hemostatic Matrix if sterile barrier package is damaged as sterility may be compromised.
- While packing a cavity for hemostasis is sometimes surgically indicated, SURGIFLO™ Hemostatic Matrix should not be used in this manner unless excess product that is not needed to maintain hemostasis is removed. When incorporated into a fibrin clot, SURGIFLO™ Hemostatic Matrix may swell up to 20% upon contact with additional fluid.
- SURGIFLO™ should not be used in conjunction with autologous blood salvage circuits.
- SURGIFLO™ should not be used in conjunction with methylmethacrylate adhesives.
- In urological procedures, SURGIFLO™ should not be left in the renal pelvis or ureters to eliminate the potential foci for calculus formation.

ADVERSE EVENTS
A total of 142 patients received SURGIFOAM® Sponge during a clinical trial comparing SURGIFOAM® Sponge to another absorbable gelatin sponge. In general, the following adverse events have been reported with the use of absorbable porcine gelatin-based hemostatic agents:

- Gelatin-based hemostatic agents may serve as a nidus for infection and abscess formation and have been reported to potentiate bacterial growth.
- Giant cell granulomas have been observed at implant sites when used in the brain.
- Compression of the brain and spinal cord resulting from the accumulation of sterile fluid have been observed.
- Multiple neurologic events were reported when absorbable gelatin-based hemostatic agents were used in laminectomy operations, including cauda equina syndrome, spinal stenosis, meningitis, arachnoiditis, headaches, paresthesias, pain, bladder and bowel dysfunction, and impotence.
- The use of absorbable gelatin-based hemostatic agents during the repair of dural defects associated with laminectomy and craniotomy operations, has been associated with fever, infection, leg paresthesias, neck and back pain, bladder and bowel incontinence, cauda equina syndrome, neurogenic bladder, impotence, and paresis.
- The use of absorbable gelatin-based hemostatic agents has been associated with paralysis, due to device migration into foramina in the bone around the spinal cord, and blindness, due to device migration in the orbit of the eye, during lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe.
- Foreign body reactions, “encapsulation” of fluid, and hematoma have been observed at implant sites.
- Excessive fibrosis and prolonged fixation of a tendon have been reported when absorbable gelatin-based sponges were used in severed tendon repair.
- Toxic shock syndrome was reported in association with the use of absorbable gelatin-based hemostats in nasal surgery.

VISTASEAL™ Fibrin Sealant (Human) - Important Safety Information
INDICATION
VISTASEAL™ is indicated as an adjunct to hemostasis for mild to moderate bleeding in adults undergoing surgery when control of bleeding by standard surgical techniques (such as suture, ligature, and cautery) is ineffective or impractical. VISTASEAL is effective in heparinized patients.

CONTRAINDICATIONS
Do not inject directly into the circulatory system.
Do not use for the treatment of severe or brisk arterial bleeding.

Do not use in patients with history of anaphylaxis or severe systemic reactions to human blood products.
Do not use VISTASEAL for spraying unless the minimum recommended distance from the applicator tip to the bleeding site can be achieved.

WARNINGS AND PRECAUTIONS
Thromboembolic events may occur if VISTASEAL is administered intravascularly. Hypersensitivity reactions can occur.
May carry a risk of transmitting infectious agents, e.g., viruses, the variant Creutzfeldt-Jakob disease (vCJD) agent and, theoretically, the Creutzfeldt-Jakob

- Fever, failure of absorption, and hearing loss have been observed when absorbable hemostatic agents were used during tympanoplasty.

Adverse Reactions to Gelatin-Based Hemostatic Agents with Thrombin
Adhesion formation and small bowel obstruction are well-known and common complications following abdominal and gynecological surgeries. Adverse events such as these, as well as inflammation and foreign body reaction including giant cell granulomas, have been reported within abdominal and gynecological surgeries in which gelatin-based and collagen-based hemostatic agents with thrombin have been used in excessive amount. Inflammation and foreign body reaction such as encapsulation of the product may mimic artifacts on radiographic images, resulting in diagnostic errors and possible reoperation. Therefore, as with other gelatin hemostatic agents, only the minimum amount of SURGIFLO™ Hemostatic Matrix needed to achieve hemostasis should be used. Once hemostasis is achieved, any excess SURGIFLO™ Hemostatic Matrix should be carefully removed.

EVITHROM™ Thrombin, Topical (Human) for Topical Use
Lyophilized Powder for Solution
EVITHROM™ is a topical thrombin indicated as an aid to hemostasis whenever oozing blood and minor bleeding from capillaries and small venules is accessible and control of bleeding by standard surgical techniques (such as suture, ligature, or cautery) is ineffective or impractical. EVITHROM™ may be used in conjunction with an Absorbable Gelatin Sponge, USP.

- Important Safety Information
- For topical use only.
 - Do not inject.
 - Apply EVITHROM™ on the surface of bleeding tissue only.
 - The amount of EVITHROM™ required depends upon the area of tissue to be treated and the method of application. In clinical studies, volumes up to 10 ml were used in conjunction with Absorbable Gelatin Sponge.
 - Do not use for the treatment of severe or brisk arterial bleeding.
 - Do not use in individuals known to have anaphylactic or severe systemic reaction to human blood products. Hypersensitivity reactions, including anaphylaxis, may occur.
 - There is a potential risk of thrombosis if absorbed systemically.
 - May carry a risk of transmitting infectious agents such as viruses and theoretically, the Creutzfeldt-Jakob disease (CJD) agent, despite manufacturing steps designed to reduce the risk of viral transmission.
 - The most common adverse reactions during clinical trial (reported in at least 2% of subjects treated with EVITHROM™) were prolonged activated partial thromboplastin time, increased INR, decreased lymphocyte count, prolonged prothrombin time, and increased neutrophil count.
 - None of the patients treated with EVITHROM™ developed antibodies to human thrombin or to human Factor V/Va. The clinical significance of these findings is unknown.

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.
US_ETH_BIOS_11036.1

disease (CJD) agent.
ADVERSE REACTIONS
The most common adverse reactions (reported in >1% of clinical trial subjects) were nausea and procedural pain.

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.

US_ETH_BIOS_110059.1

EVARREST® Fibrin Sealant Patch - Important Safety Information
INDICATIONS AND USAGE
EVARREST® is a fibrin sealant patch indicated for use with manual compression as an adjunct to hemostasis in adult patients undergoing surgery, when control of bleeding by conventional surgical techniques (such as suture, ligature, and cautery) is ineffective or impractical.

LIMITATIONS FOR USE

- Cannot be used in place of sutures or other forms of mechanical ligation in the treatment of major arterial or venous bleeding.
- Not for use in children under one month of age
- Laparoscopic and other minimally invasive surgeries where manual compression would be difficult to achieve.

DOSAGE AND ADMINISTRATION
For topical use only

- Determine the number of patches to be applied based upon the surface area and anatomic location of the bleeding tissue to be treated.
- Keep the patch dry until use.
- Place the powdery (active) side of the patch on the surface of tissue.
- Apply immediate manual compression over the entire surface of the patch and maintain contact pressure for 3 minutes to control the bleeding.

DOSAGE FORMS AND STRENGTHS
EVARREST® Fibrin Sealant Patch consists of human fibrinogen and human thrombin embedded in a flexible composite patch component. The active side is powdery, and the non-active side has an embossed wave pattern.
Each 2 x 4 inch (5.1 x 10.2 cm) absorbable patch contains:

- 55.5 mg per square inch (8.6 mg per square cm) human fibrinogen
- 241.9 Units per square inch (37.5 Units per square cm) human thrombin

CONTRAINDICATIONS

- Do not use to treat bleeding from large defects in arteries or veins.
- Do not apply intravascularly.
- Do not use in individuals known to have anaphylactic or severe systemic reaction to human blood products.

WARNINGS AND PRECAUTIONS

- Thrombosis can occur if absorbed systemically. Apply topically to the bleeding site only.
- Can cause hypersensitivity reactions including anaphylaxis.
- Avoid application to contaminated areas of the body or in the presence of active infection. Infection can occur.
- EVARREST® contains oxidized regenerated cellulose which adheres to bleeding surfaces. Inadvertent adhesions can occur.
- Avoid use in, around, or in proximity to, foramina in bone or areas of bony confine where swelling may cause compression.
- Use the least number of patches required to cover the entire bleeding area.
- May carry a risk of transmitting infectious agents, e.g., viruses, the variant Creutzfeldt-Jakob disease (vCJD) agent and, theoretically, the Creutzfeldt-Jakob disease (CJD) agent.

***The bleeding situations identified reflect customer insights/market research on optimal adjunctive hemostat utilization. The product solutions should only be used in accordance with their instructions for use. Product recommendations should not supplant medical judgment. Surgeon preference, experience, and patient needs may dictate alternate technique. Review all relevant precautions, especially the indications, contraindications, warnings, and information for use. Please see package inserts for Full Prescribing Information. The visual does not reflect any sequential order in use.**

References: 1. Hemostasis Optimization Program and Next Steps Meeting, dated 06/30/2014.